

50_{ml} Centrifuge Tubes

15_{ml} Tubes





- ✓ Caps easy to close with one hand
- ✓ Autoclavable

- ✓ Autoclavable
- ✓ Caps easy to close with one hand
- ✓ Available in e-beam pre-sterilized format
- ✓ Printed graduations and writing area for sample identification
- ✓ Made from Polypropylene Copolymer for superior strength

Watch our video review



Part	Capacity	Description	Tube Dimensions H x W	Cap Dimensions H x W	Material
 WE4G	15ml	Centrifuge Tube, Grad. Conical, Blue Cap, Pre-sterilized	119 x 18mm	11 x 20mm	P/P
Packaging 500 bulk					
 WE2G	50ml	Centrifuge Tube, Grad. Conical, Blue Cap, Pre-sterilized	116 x 30mm	13 x 33mm	P/P
Packaging 500 bulk					

For more information or to request samples, please contact Pacific Laboratory Products at: sales@pacificlab.com.au

Distributed by



Australian Head Office
PO BOX 465
Blackburn, VIC 3130

Free Call: 1800 723 405
Ph: (03) 9845 0300
Fax: + 61 3 9845 0350
e-mail: sales@pacificlab.com.au

Connect with us:

