

Introducing



Gel Documentation Systems & Horizontal Gel Boxes from Axygen®

Ideal for advanced genomic workflow

Pacific Laboratory Products introduces the latest Gel Documentation Systems from US manufacturer Axygen®. With various models available to select according to needs - from basic to advanced, the system allows for visualization of DNA and protein stained within a gel photo documentation in a laboratory setting.


Key advantages of the system:

- One piece integrated digital darkroom
- 5.4 megapixel high resolution camera
- Built-on windows tablet PC with 4 click data acquisition
- Blue light LED for imaging of DNA dyes at 470nm
- 30,000 hour dual wavelength UV transilluminator
- Plug and play devices
- Easy to use



Learn more! Watch the Video

<http://bit.ly/2oABXbQ>



With a small footprint these units are presented with an integrated cabinet "self contained darkroom" of 28lbs, 22 x 12 x 15in.

Warranty

The system is available with a two years warranty including the UV bulbs, and are made in the USA.



The system allows users to obtain an image in just four clicks:
Select darkroom ▶ light source ▶ auto expose ▶ capture image... it is that simple!

[Learn more! Click here](#)



The Axygen® Horizontal Gel Boxes are ideal to work with reagents to allow the separation of DNA based on size in an Agarose matrix. These are then imaged on the Gel Doc Systems.

[Learn more! Click here](#)

There are four sizes available for these horizontal tanks.

Each comes with two size gel trays that allow for external gel casting.

The three large sizes come with recirculation ports for use as needed.



The Axygen Gel Boxes are an affordable alternative and a great match when used with the Axygen Gel Doc Systems. Axygen Gel Boxes speed up workflow by allowing work to be done in parallel instead of serially.

For more information simply contact one of our product experts at sales@pacificlub.com.au

Distributed by Pacific Laboratory Products | PO BOX 465 Blackburn, VIC 3130 | Free Call: 1800 723 405 Ph: (03) 9845 0300 | e-mail: sales@pacificlub.com.au

www.pacificlub.com.au

[f /PacifLabP](#) [@PacifLabP](#)