

Biopette E™ Electronic Pipette

The Biopette E is a multifunctional pipette that is fully motorized, eliminating hand strain and allowing efficient multidispensing.

A high precision stepping motor controls piston action, resulting in exceptional accuracy and precision and minimal pipetting forces. Accuracy is further enhanced by a state-of-the-art microprocessor which controls the action of the stepping motor. The calibration of the unit is checked each time an operation is performed.

The Biopette E features an adjustable tip ejector along with a universal tip shaft designed for compatibility with a wide variety of manufacturers' tips. The lower portion of the pipette is easily removed for autoclaving.

Each pipette is individually tested and supplied with a certificate of quality. Also included is a shelf clip, Lithium ion battery and battery charger. An optional charger set is available for charging the battery outside of the pipette.

- Comfortable and effortless pipetting
- Efficiency at an affordable price
- Exceptional accuracy and precision
- Fully motorized piston drive
- User-friendly operator interface



Special Pricing

Last Unit in Stock!

\$395.00*



SPECIFICATIONS

Volume Range	Tip	Increments	Accuracy	Precision
0.5-10 µl	10 µl	0.1 µl	±4.0 to ±1.0%	<2.5 to <0.4%
2-20 µl	10 µl*	0.1 µl	±5.0 to ±1.0%	<2.0 to <0.3%
10-200 µl	200 µl	1.0 µl	±2.0 to ±0.6%	<1.0 to <0.15%
100-1,000 µl	1,000 µl	1.0 µl	±1.5 to ±0.5%	<0.6 to <0.15%

*Most manufacturers' 10 µl tips will accurately accommodate 20 µl.

Catalog No.

Description

P3608-200



BioPette Digital, 10 - 200 µL, 8 - Channel

* Price quoted in AUD exclude GST and delivery charges. Offer excludes spare parts, servicing, freight, shipping and handling and is available while stock last Only. Images displayed on this publication are a real representation of the product on sale.

Distributed by



www.pacificlab.com.au

Australian Head Office

PO BOX 465
Blackburn, VIC 3130

Connect with us:

Free Call: 1800 723 405

Ph: (03) 9845 0300

Fax: + 61 3 9845 0350

e-mail: sales@pacificlab.com.au



/PacificLabP



@PacificLabP