

NEW ERMA TRIMMING HANDLE

for Microtome Blades

Part No. 08-660-2

[Shop](#)

**Introductory
Offer!**




Unique Snap-Lock Mechanism for safety

- This new safety designed trimming holder can be used with all brands of disposable microtome blades to cut-up.
- The Trimming Holder comes with a unique built-in Snap-Lock mechanism which holds the blade firmly in place.
- The handle is heat resistant, chemically resistant and can be autoclaved.

Mount the blade in just few simple steps:

Press to release and open the trimming holder, mount the blade, then close to secure the blade firmly.



Part No.	Description	Qty	Price AUD Exc. GST
08-660-2	Trimming Handle For Microtome Blade	5 per pack	Introductory Offer! \$ 95.00 Shop



The trimming holder is fully compatible and can be used with the ERMA Microtome Blades.
(Microtomes blades sold separately)

ERMA INC.
TOKYO JAPAN

Distributed by



www.pacificlab.com.au

Australian Head Office

PO BOX 465
Blackburn, VIC 3130

Connect with us:

Free Call: 1800 723 405

Ph: (03) 9845 0300

Fax: + 61 3 9845 0350

e-mail: sales@pacificlab.com.au

[f](#) /PacificLabP [t](#) @PacificLabP