

GENESIG

NEW Coronavirus Test kit

Now available at Pacific Laboratory Products

COVID-19 for specific detection of 2019-nCoV outbreak Test kit excludes bat SARS-related coronavirus

Primerdesign are the first to design a target specific assay for the surging outbreak of Coronavirus COVID-19. This strain is highly contagious, known to cause fatalities and is spreading from person to person, causing global impact.

KEY BENEFITS

- Clinically Diagnostic Test (CE-IVD)
- Genesig Real time PCR COVID-19 Kit (CE)
- 100% detection profile
- Rapid results within 2 hours
- Specific to the latest Coronavirus outbreak
- Room temperature – fast shipping



Genesig nCoV-2019 test kit has been tested and is compatible with the following real-time qPCR instruments:
Roche Light cycler 480, Bio-Rad CFX 96, Applied Biosystems 7500.

Patient samples: Oropharyngeal Swab, nasopharyngeal swab, sputum.

Genesig™ provides the world's broadest menu of qPCR detection kits, fast development of new assays on demand. Additionally, we design, validate and manufacture qPCR kits, master mixes, variety of qPCR controls, lyophilised reagents, and qPCR instruments.

Pacific Laboratory Products
sales@pacificlab.com.au
Free Call: 1800 723 405
Ph: 03 9845 0300
www.pacificlab.com.au
www.genesig.com

/PacificLabP @PacificLabP