

**NEW**

- Powder free ✓
- Long cuff ✓
- Fingertip textured ✓
- Lab Chem tested ✓
- Non-sterile ✓
- Ambidextrous ✓



Samples available on request

# gloveon Nitrile Examination Gloves

Developed for maximum comfort

NON STERILE CE  
Extended Cuff



## MATADOR 8mil Nitrile Examination Gloves

### Superior Protection



50 Gloves Per box


- Developed to be used in harsh environments
- Ideal for industrial applications where extra thickness is required
- Nitrile exclusive formula offers unmatched elasticity and low tearing rate



100 Gloves Per box

## VIGOR

- Robustness and hardness provides extra protection well over the wrist
- Unique high visibility colour, suitable for industrial environments
- Clinically tested to protect against a wide range of chemicals complying with various safety standards



- ☑ Nitrile
- ☑ Eureka
- ☑ Celeste
- ☑ Celeste LC
- ☑ Coats
- ☑ Maverick
- ☑ Maverick LC
- ☑ Matador
- ☑ Paladin
- ☑ Vigor
- ☑ Nitrile Sterile
- ☑ Aegis (Sterile)
- ☑ Latex
- ☑ Ridley
- ☑ Dispensers

[Learn more!](#)  
[Watch the video](#)

Distributed By



Free Call: 1800 723 405  
Ph: (03) 9845 0300  
email: [sales@pacificlab.com.au](mailto:sales@pacificlab.com.au)

[www.pacificlab.com.au](http://www.pacificlab.com.au)

### Find the right glove for you!

Check out the full range, see the chart on next page


## Nitrile Latex Vinyl




11 different types of gloves to select from

[Catalogue](#)

## Range

**Eureka** 



**Celeste** 



**Celeste LC** 



**Coats** 



**Maverick** 




**Maverick LC** 




UL CERTIFIED

**Paladin** 





**Matador** 



**Vigor** 



**Aegis** 



STERILE - R

**Ridley** 



## Dispensers

Range	Type	Textured	Powder free	Cuff	Length (mm)	Thickness (mm)	Sterile	UL Certified	Ambi-dextrous	Natural Rubber Made	Colour	Applications	Ordering Code	Qty Per Box	Boxes x Ctn
<b>Eureka</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Standard Cuff	≥ 230	Palm ≥ 0.05 Finger ≥ 0.07	Non-Sterile	✗	✓	✗	Dark Blue	• General laboratory	NTR61XS NTR61SS NTR61MM NTR61LL NTR61XL	300 300 300 300 270	10
<b>Celeste</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Standard Cuff	≥ 230	Palm 0.07 ± 0.02 Finger 0.10 ± 0.02	Non-Sterile	✗	✓	✗	Violet Blue	• Delicate procedures • High tactile sensitivity requirements • Chemical protection	NTR58XS NTR58SS NTR58MM NTR58LL NTR58XL	200 200 200 200 180	10
<b>Celeste LC</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Long Cuff	≥ 280	Palm 0.08 ± 0.02 Finger 0.12 ± 0.02	Non-Sterile	✗	✓	✗	Violet Blue	• Delicate procedures • High tactile sensitivity requirements • Chemical protection	NTR56SS NTR56MM NTR56LL NTR56XL	150	10
<b>Coats</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Standard Cuff	≥ 240	Palm 0.07 ± 0.02 Finger 0.09 ± 0.02	Non-Sterile	✗	✓	✗	Light Grey	• General laboratory • Special for users with sensitive skin	CTS38XS CTS38SS CTS38MM CTS38LL CTS38XL	200 200 200 200 180	10
<b>Maverick</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Standard Cuff	≥ 230	Palm ≥ 0.08 Finger ≥ 0.11	Non-Sterile	✗	✓	✗	Aqua Blue	• Chemistry • Harsh environments • Industrial	NTR28XS NTR28SS NTR28MM NTR28LL NTR28XL	200 200 200 200 180	10
<b>Maverick LC</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Long Cuff	≥ 280	Palm 0.09 ± 0.02 Finger 0.14 ± 0.02	Non-Sterile	✓	✓	✗	Aqua Blue	• Chemistry • Harsh environments • Industrial	NTR26XS NTR26SS NTR26MM NTR26LL NTR26XL	100	10
<b>Paladin</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Long Cuff	≥ 280	Palm 0.10 ± 0.02 Finger 0.13 ± 0.02	Non-Sterile	✗	✓	✗	Cobalt Blue	• Chemistry • Corrosive environments • Industrial	NTR36XS NTR36SS NTR36MM NTR36LL NTR36XL NTR36XXL NTR36XXXL	100	10
<b>Matador</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Long Cuff	≥ 280	Palm 0.20 ± 0.02 Finger 0.24 ± 0.02	Non-Sterile	✗	✓	✗	Sky Blue	• Harsh environments • Industrial • Where extra thickness is required	NTR18XS NTR18SS NTR18MM NTR18LL NTR18XL	50	10
<b>Vigor</b> 	<b>Nitrile</b>	Fingertip Textured	✓	Long Cuff	≥ 280	Palm 0.13 ± 0.02 Finger 0.17 ± 0.02	Non-Sterile	✓	✓	✗	Orange	• Harsh environments • Industrial • Where extra protection is required	NTR27SS NTR27MM NTR27LL NTR27XL	100	10
<b>Aegis</b> 	<b>Sterile Nitrile</b>	Fingertip Textured	✓	Long Cuff	≥ 280	Palm 0.08 ± 0.02 Finger 0.12 ± 0.02	Sterile-R	✗	✓	✗	Violet Blue	• Animal handling • Clinicians • Special laboratory	SNT56SS SNT56MM SNT56LL SNT56XL	50 Pairs	4
<b>Ridley</b> 	<b>Latex</b>	Fully Textured	✓	Standard Cuff	≥ 230	Palm 0.13 ± 0.03 Finger 0.14 ± 0.03	Non-Sterile	✗	✓	✓	Natural Colour	• General laboratory • Life sciences	CHL11XS CHL11SS CHL11MM CHL11LL CHL11XL	100	10
<b>Dispensers</b>	<b>Acrylic</b>		-			94D x 139L x 225W (mm) 100D x 405L x 230W (mm) 100D x 533L x 230W (mm) 150D x 200W (mm) 150D x 200W (mm)			-		Clear	• Holds one box • Holds three boxes • Holds four boxes • Multipurpose • Multipurpose	DS1 DS30 DS40 DSM10 DSM20	each	-