

LABOCAP



GERBER
INSTRUMENTS
Since 1892 



LABOCAP with handle

- Material:** Aluminium, colored anodized.
Caps close vessels sterile, but not hermetically.
The Springs clamp to the outer vessel wall.
- Advantages:** rust proof, washable*, reusable, heat proof, sterilizable, various colors for the recognition, very good adhesion, airflow is guaranteed.
(green, yellow, violet and black are special colors with supplement)

For glasses with Ø mm	Pc.	Blue	Red	Silver	Green	Yellow	Violet	Black
15/16	100	9.231403	9.231413	9.231345	9.231327	9.231326	9.231329	9.231328
17/18	100	23.21043	23.21042	9.231347	9.231331	9.231330	9.231333	9.231332



LABOCAP without handle

- Material:** Aluminium, colored anodized.
Caps close vessels sterile, but not hermetically.
The Springs clamp to the outer vessel wall.
- Advantages:** rust proof, washable*, reusable, heat proof, sterilizable, various colors for the recognition, very good adhesion, airflow is guaranteed.
(green, yellow, violet and black are special colors with supplement)

For glasses with Ø mm	Pc.	Blue	Red	Silver	Green	Yellow	Violet	Black
9/10	100	9.231421	9.231441	9.231369	23.21210	23.21212	23.21213	23.21214
12/13	100	9.231422	9.231442	9.231372	9.231463	9.231460	9.231335	9.231334
14/15	100	9.231423	9.231443	9.231374	23.21215	23.21216	23.21217	23.21218
15/16	100	9.231424	9.231444	9.231375	9.231464	9.231461	9.231336	9.231466
17/18	100	9.231425	9.231445	9.231377	9.231465	9.231462	9.231338	9.231337
19/20	100	9.231426	9.231446	9.231379	9.231340	9.231339	9.231342	9.231341
21/23	100	9.231427	9.231447	9.231381	23.21219	23.21220	23.21221	23.21222
24/26	100	9.231428	9.231448	9.231385	23.21223	23.21224	23.21225	23.21226
28/30	100	9.231429	9.231449	9.231388	23.21227	23.21228	23.21229	23.21230



ERLENMEYERCAP without handle

- Material:** Aluminium, colored anodized.
Caps for culture flasks (without beaded rim) and nutrient media bottles with straight neck (Ø 38 mm). The Springs clamp to the outer vessel wall.
- Advantages:** rust proof, washable*, reusable, heat proof, sterilizable, various colors for the recognition, very good adhesion, airflow is guaranteed.
(green, yellow, violet and black are special colors with supplement)

For glasses with Ø mm	Pc.	Blue	Red	Silver	Green	Yellow	Violet	Black
37/39	10	9.231457	9.231458	9.231456	23.21231	23.21232	23.21233	23.21234

*Cleaning recommendations: The colored anodized coating on the aluminum cap is not resistant to alkaline cleaners. For cleaning will a mild neutral detergent recommended (detergent).

GERBER
INSTRUMENTS
Since 1892

Distributed by
 PLP
PACIFIC LABORATORY PRODUCTS

email: sales@pacificlab.com.au
Free Call: 1800 723 405
Ph: (03) 9845 0300
www.pacificlab.com.au