



Frontier 5515 Micro Centrifuge

Special Price
\$4,360.00*

12
Months
Warranty

Combo Deal

** For limited time only!*

5515 Microliter Benchtop Centrifuge

Price Drop!

High-Speed Microliter Centrifuge for Liquid Separation Applications

Compatible with a variety of tube volumes and capacities, the 5515 Benchtop centrifuge provides high-speed operation for spin-down applications.

- Maximum g-force of 21,953 RCF
- Powerful induction motor reaches speeds up to 15,200 RPM in less than 16 seconds
- 10 acceleration/deceleration settings



ID: 30130866



ID: 30130884



ID: 30130871

[find out more](#)

Combo Includes Centrifuge, Rotor & Adapters

Combo Deal Includes:

30130866 / FC5515 Centrifuge + 30130871 / 24x1.5/2ml sealable rotor + 30130884 / Adapters for 0.2/0.4 Tubes (4 x packets of 6=24)



An intuitive control panel with a rotary knob allows the centrifuge to be operated easily.

A waterproof keypad provides protection that helps significantly increase the lifespan of the unit.



A powerful brushless induction motor, 10 acceleration and deceleration settings, and memory for up to 99 commonly-used parameters are only a few features of this unit.

The Centrifuge is compatible with seven different rotors (additional rotors are optional).

Please quote Promotion Code: **OHA-1218-1** when ordering by email or telephone. Combo deal not available online please call/email.

Terms and Conditions

* Prices quoted in AUD exclude GST and delivery charges.

* Offer is valid until the 31st of Dec 2018 Only.

* Promotion cannot be used in conjunction with any other offer. Offer exclude spare parts, servicing, freight, shipping and handling and is available while stocks last Only, no rain-checks. Offer is not available for re-sellers. Prices apply to stock "On Hand" Only. Pacific Laboratory Products Pty Ltd reserves the right to modify or withdraw this offer at any time.

Images displayed on this publication are a graphic representation of the product only and may vary from the item itself.

www.pacificlab.com.au

Australian Head Office

PO BOX 465
Blackburn, VIC 3130
Free Call: 1800 723 405
Ph: (03) 9845 0300
Fax: + 61 3 9845 0350

e-mail: sales@pacificlab.com.au

[f](#) /PacificLabP [t](#) @PacificLabP