



# Benchtop High-Speed MicroCentrifuge

## TG1650-WS



### Introducing our **NEW** Micro-Centrifuge alternative

An affordable and reliable unit ideal for molecular research applications



### Introductory Offer

- Micro-processor controlled
- Brushless induction drive motor
- Reliable, stable and low noise operation with high speed accuracy
- Includes self-locking device, over-speed safety protection, and automatic alarm
- The real-time conversion of Speed and RCF is convenient for operation
- Programmable operation, data can be set according to demands and stored automatically
- Intuitive control panel. User-friendly interface with LED or LCD display for simple and convenient operation
- 9 acceleration and deceleration selectable profiles

**\$1,650.00\***

Bundle: with 1.5/2ml x 24 Rotor

ID: TG1650-WS

**12**  
Months  
Warranty

Specifications			
Model	TG1650-WS	Time Setting	1min ~ 99min
Max Speed	16,500 r/min	Noise	<65dB(A)
Max RCF	23,669 x g	Power supply	220V/50Hz 10A
Max Capacity	100ml x 4	Dimension (LxWxH)	330mm x 420mm x 280mm
Speed Accuracy	± 30r/min	Packing size (LxWxH)	430mm x 520mm x 390mm
		Net weight	20kg

### Included Rotor



**TGL-17M-R3**  
Capacity: 1.5 / 2ml x 24  
Speed: 13,200 r/min  
Centrifugal force: 16,110 x g  
Inc. Set Adapters: 0.2ml, 0.5ml

### Optional Compatible Rotors



**TGL-17M-R1**  
Capacity: 12 x 1.5 / 2ml  
Speed: 16,500 r/min  
Centrifugal force: 18,930 x g  
Inc. Set Adapters: 0.2ml, 0.5ml



**TGL-17M-R2**  
Capacity: 12 x 5ml  
Speed: 13,000 r/min  
Centrifugal force: 11,400 x g  
Inc. Set Adapters: 0.5ml, 1.5ml



**TGL-17M-R5**  
Capacity: 48 x 0.5ml  
Speed: 16,000 r/min  
Centrifugal force: 23,669 x g



**TGL-17M-R6**  
Capacity: 3 x (8 x 0.2ml PCR)  
Speed: 15,000 r/min  
Centrifugal force: 13,800 x g

\* Price quoted in AUD excludes GST and delivery charges.

Distributed by



[www.pacificlab.com.au](http://www.pacificlab.com.au)

**Australian Head Office**

PO BOX 465  
Blackburn, VIC 3130

Connect with us:

/PacificLabP @PacificLabP

Free Call: 1800 723 405

Ph: (03) 9845 0300

Fax: + 61 3 9845 0350

e-mail: [sales@pacificlab.com.au](mailto:sales@pacificlab.com.au)